

# Emerald™ Design Fitting Certification Course

## *Test Booklet*



Let's start the Emerald Design Certification process. Here are the directions for completing your **Answer Sheet** (last page).

- Be sure to completely fill-out the personal information on the answer sheet & forward to ACL. PLEASE PRINT LEGIBLY
- Choose the **best** answer from those offered.
- Mark your choice by circling your answer.
- If you decide to change your answer draw an "X" through the answer you do not want and circle your next choice.

### Question 1

In which direction are the fluids compressed or redistributed when a ortho-k lens is worn?

- a. From the centre toward the periphery
- b. From the anterior toward the posterior
- c. From the periphery toward the centre
- d. From the posterior toward the anterior

### Question 2

What is the value of topography in the course of ortho-k treatment?

- a. It has limited or no value in the ortho-k treatment process.
- b. It is only required for choosing the initial ortho-k shaping trial lens in some fitting systems
- c. Whether used for choosing the initial ortho-k shaping lens or not, topography is necessary to monitor the entire ortho-k treatment effect process
- d. It is used in conjunction with keratometry to determine the initial ortho-k shaping trial lens or the first shaping lens to be ordered.



### Question 3

**What is the function of the second fitting/reverse curve of the ortho-k shaping lens?**

- a. To provide the optical correction
- b. To promote movement of the ortho-k shaping lens across the cornea
- c. To aid the forces controlling ortho-k shaping lens centration
- d. To allow fluid and cell “movement” which promotes mid-peripheral corneal thickening

### Question 4

**What is the function of the third or alignment curve of the ortho-k shaping lens?**

- a. To promote edge lift
- b. To stabilize the ortho-k shaping lens and promote optimal centration
- c. To link the back optic zone radius (BOZR) to the reverse curve
- d. To minimize movement of the ortho-k shaping lens

### Question 5

**What is the approved maximum “against-the-rule” astigmatism that can be attempted with Emerald ortho-k fitting?**

- a. -0.75 Diopters
- b. -2.75 Diopters
- c. -0.25 Diopters
- d. -1.75 Diopters

### Question 6

**What is the approved amount of myopic reduction that can be achieved with Emerald ortho-k treatment?**

- a. 10.00 Diopters
- b. 8.00 Diopters
- c. 5.00 Diopters
- d. 2.00 Diopters



### Question 7

**Why may people with a large pupil size experience problems with ortho-k?**

- a. The treatment zone will not be large enough and will cause visual problems
- b. The treatment zone will be too large and will cause myopic regression
- c. The treatment zone will be too large and will cause night-time blur
- d. The treatment zone will not be large enough and will cause lens decentration

### Question 8

**Why is significant internal or lenticular astigmatism a problem with ortho-k fitting?**

- a. Ortho-k shaping lenses are only produced with a back-surface cylinder correction
- b. Ortho-k always causes a decrease in the cylinder correction
- c. Ortho-k shaping lenses only impact corneal astigmatic changes
- d. Ortho-k shaping lenses are difficult to manufacture

### Question 9

**Which patient symptoms with vision shaping treatment are of most concern to the contact lens practitioner?**

- a. Comfortable wear of the ortho-k shaping lens and facial redness
- b. Lid twitching and loss of lashes
- c. Pain, lacrimation, redness, photophobia
- d. Inability to sleep and sweating

### Question 10

**What is the most important common observation that the wearer should make prior to attempting removal the ortho-k shaping lenses after waking?**

- a. Check for mucus build-up in the nasal canthus
- b. Check for adherence of the ortho-k shaping lens to the cornea
- c. Check for redness around the limbus
- d. Check for excessive discomfort in the eye



### Question 11

Which of the following patients would be considered a *poor candidate* to be fitted with the Emerald Design?

- a. -3.00 -1.00 X 180
- b. -0.50 -1.75 X 180
- c. -5.00 -0.50 X 90
- d. -2.50 -1.50 X 180

### Question 12

The recommended data needed for designing an Emerald lens is?

- a. Refraction, accurate "K's", and VD
- b. Refraction, accurate "K's", and HVID
- c. Refraction, accurate "K's", and IOP
- d. Refraction, accurate "K's", and VST

### Question 13

Which description best identifies the Emerald Ortho-k Design?

- a. A four-zone reverse geometry lens available in three diameters
- b. An ortho-k design with spherical central zone and toric periphery
- c. Aspheric back surface with a front surface add
- d. A bi-aspheric lens with two variable parameters

### Question 14

The initial recommended Emerald lens diameter is?

- a. 11.5 mm
- b. 10.2 mm
- c. 10.0 mm
- d. 10.6 mm



### Question 15

The Emerald design typically uses an Alignment Curve radius that is?

- a. Steeper than the Flat "K" radius
- b. Flatter than the Flat "K" radius
- c. Equal to the Flat "K" radius
- d. Aspheric

### Question 16

What parameter change would likely be required if an Emerald Design lens decentres horizontally?

- a. Increase the overall lens diameter
- b. Decrease the overall lens diameter
- c. Steepen the Base (BOZR) curve
- d. Steepen the Peripheral Curve

### Question 17

A Subtractive or Difference map is used to determine?

- a. The optimum design for the back surface of the ortho-k shaping lens
- b. The difference between the pre and post-fit cornea
- c. The sagittal height of the cornea
- d. The baseline eccentricity value of the cornea

### Question 18

A "Central Island" topography map would indicate that the Emerald lens?

- a. Has the proper sagittal height
- b. Has a sagittal height that is too low
- c. Has a sagittal height that is too high
- d. Needs to have the diameter changed



### **Question 19**

**A “Smiley Face” topography map would indicate that the Emerald lens?**

- a. Is centering ideally
- b. Is centering high on the cornea during sleep
- c. Is centering low on the cornea during sleep
- d. Is decentred nasally on the cornea

### **Question 20**

**A “Bull’s Eye” topography map would indicate that the Emerald lens?**

- a. Has the proper sagittal height
- b. Has a sagittal height that is too low
- c. Has a sagittal height that is too high
- d. Needs to have the diameter changed

***Thank You for Completing the Test***



## CERTIFICATION TEST ANSWER SHEET

Date:

First Name:

Last Name:

Practice Name:

Address:

State:

Postcode:

Phone:

E-mail address:

**Please circle one answer for each question.**

(If you decide to change an answer, please cross out the incorrect answer with an "X" and circle new choice.)

- |     |   |   |   |   |     |   |   |   |   |
|-----|---|---|---|---|-----|---|---|---|---|
| 1.  | A | B | C | D | 11. | A | B | C | D |
| 2.  | A | B | C | D | 12. | A | B | C | D |
| 3.  | A | B | C | D | 13. | A | B | C | D |
| 4.  | A | B | C | D | 14. | A | B | C | D |
| 5.  | A | B | C | D | 15. | A | B | C | D |
| 6.  | A | B | C | D | 16. | A | B | C | D |
| 7.  | A | B | C | D | 17. | A | B | C | D |
| 8.  | A | B | C | D | 18. | A | B | C | D |
| 9.  | A | B | C | D | 19. | A | B | C | D |
| 10. | A | B | C | D | 20. | A | B | C | D |

Please submit completed answer sheet to Australian Contact Lenses via email [aclorders@gelflex.com](mailto:aclorders@gelflex.com) or Fax 03 9793 1635.