

## Extended Diopter Conversion Chart

DIOPTER	RADIUS	DIOPTER	RADIUS	DIOPTER	RADIUS	DIOPTER	RADIUS	DIOPTER	RADIUS
23.00	14.67	35.00	9.64	41.00	8.23	47.00	7.18	53.00	6.36
24.00	14.06	35.12	9.61	41.12	8.20	47.12	7.16	53.12	6.35
25.00	13.50	35.25	9.57	41.25	8.18	47.25	7.14	53.25	6.33
26.00	12.98	35.37	9.54	41.37	8.15	47.37	7.12	53.37	6.32
27.00	12.50	35.50	9.50	41.50	8.13	47.50	7.10	53.50	6.30
27.50	12.27	35.62	9.47	41.62	8.10	47.62	7.08	53.62	6.29
28.00	12.05	35.75	9.44	41.75	8.08	47.75	7.06	53.75	6.27
29.00	11.63	35.87	9.40	41.87	8.06	47.87	7.05	53.87	6.26
30.00	11.25	36.00	9.37	42.00	8.03	48.00	7.03	54.00	6.25
30.12	11.20	36.12	9.34	42.12	8.01	48.12	7.01	54.12	6.23
30.25	11.15	36.25	9.31	42.25	7.98	48.25	6.99	54.25	6.22
30.37	11.11	36.37	9.28	42.37	7.96	48.37	6.97	54.37	6.20
30.50	11.06	36.50	9.24	42.50	7.94	48.50	6.95	54.50	6.19
30.62	11.02	36.62	9.21	42.62	7.91	48.62	6.94	54.62	6.17
30.75	10.97	36.75	9.18	42.75	7.89	48.75	6.92	54.75	6.16
30.87	10.93	36.87	9.15	42.87	7.87	48.87	6.90	54.87	6.15
31.00	10.88	37.00	9.12	43.00	7.84	49.00	6.88	55.00	6.13
31.12	10.84	37.12	9.09	43.12	7.82	49.12	6.87	55.12	6.12
31.25	10.80	37.25	9.06	43.25	7.80	49.25	6.85	55.25	6.10
31.37	10.75	37.37	9.03	43.37	7.78	49.37	6.83	55.37	6.09
31.50	10.71	37.50	9.00	43.50	7.75	49.50	6.81	55.50	6.08
31.62	10.67	37.62	8.97	43.62	7.73	49.62	6.80	55.62	6.06
31.75	10.63	37.75	8.94	43.75	7.71	49.75	6.78	55.75	6.05
31.87	10.59	37.87	8.91	43.87	7.69	49.87	6.76	55.87	6.04
32.00	10.54	38.00	8.88	44.00	7.67	50.00	6.75	56.00	6.03
32.12	10.50	38.12	8.85	44.12	7.65	50.12	6.73	56.25	6.00
32.25	10.46	38.25	8.82	44.25	7.62	50.25	6.71	56.50	5.97
32.37	10.42	38.37	8.79	44.37	7.60	50.37	6.70	56.75	5.95
32.50	10.38	38.50	8.76	44.50	7.58	50.50	6.68	57.00	5.93
32.62	10.34	38.62	8.73	44.62	7.56	50.62	6.66	57.25	5.90
32.75	10.30	38.75	8.71	44.75	7.54	50.75	6.65	57.50	5.88
32.87	10.26	38.87	8.68	44.87	7.52	50.87	6.63	57.75	5.85
33.00	10.22	39.00	8.65	45.00	7.50	51.00	6.61	58.00	5.83
33.12	10.19	39.12	8.62	45.12	7.48	51.12	6.60	58.25	5.80
33.25	10.15	39.25	8.59	45.25	7.45	51.25	6.58	58.50	5.77
33.37	10.11	39.37	8.57	45.37	7.43	51.37	6.57	58.75	5.75
33.50	10.07	39.50	8.54	45.50	7.41	51.50	6.55	59.00	5.72
33.62	10.03	39.62	8.51	45.62	7.39	51.62	6.53	59.25	5.69
33.75	10.00	39.75	8.49	45.75	7.37	51.75	6.52	59.50	5.67
33.87	9.96	39.87	8.46	45.87	7.35	51.87	6.50	59.75	5.65
34.00	9.92	40.00	8.43	46.00	7.33	52.00	6.49	60.00	5.63
34.12	9.89	40.12	8.41	46.12	7.31	52.12	6.47	60.25	5.60
34.25	9.85	40.25	8.38	46.25	7.29	52.25	6.45	60.50	5.58
34.37	9.82	40.37	8.36	46.37	7.27	52.37	6.44	60.75	5.56
34.50	9.78	40.50	8.33	46.50	7.25	52.50	6.42	61.00	5.53
34.62	9.74	40.62	8.30	46.62	7.23	52.62	6.41	61.25	5.51
34.75	9.71	40.75	8.28	46.75	7.21	52.75	6.39	61.50	5.49
34.87	9.67	40.87	8.26	46.87	7.20	52.87	6.38	61.75	5.47

<b>Diopter to Radius Conversion</b> $\frac{337.5}{D} = \text{mm}$	<b>Radius to Diopter Conversion</b> $\frac{337.5}{\text{mm}} = D$
--	--

